Field Safety Notice

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Dear Customers,
Chers clients

Vital Healthcare has identified the possibility of abnormal device behaviors of the NovaLine® Tubing Sets for Hemodialysis listed above.

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Under certain rare circumstances, the device may exhibit the following behaviors:

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Vital Healthcare has taken the customer feedback very seriously and conducted a thorough and detailed investigation to better understand the circumstances of the reported events. Risk analysis executed based on the relevant risk management files has underlined that residual risks are acceptable after taking necessary corrective and preventive actions, minimizing the occurrence of potential risks. Considering that no direct health risk / threat against patients was identified and no serious health threat happened to involved patients, no product recall is needed at present.

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The **Annex 1** provides the “Action to be Taken by Customers / Users” on how to identify affected devices and instructions on actions to be taken. To address these actions, Vital Healthcare is providing an Instructions for Use (IFU) Addendum in **Annex 2**. This Addendum describes how users can check the suitability of the device for the intended use, providing measures to prevent, or methods to allow the user to correct the abnormal devices.

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Vital Healthcare apologizes for any inconvenience this may cause you. If you have questions regarding this notification or need any further information or support, please contact Baxter Quality Department on +33 (0) 1 34 61 54 93.

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Vital Healthcare Sdn. Bhd. confirms that the competent authority has been informed of this Field Safety Notice.

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Transmission of this Field Safety Notice:

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Please send this notice to all concerned persons of your organization and / or organization where the products were transferred.

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Please check the application of this notice to ensure the effectiveness of corrective measures recommended.

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Cordialement,

Vital Healthcare Sdn. Bhd.

Annex 1: Action to be Taken by Customers / Users

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Annex 2: Additional Instructions for Use

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**Action to be Taken by
Customers / Users**

Operator must check the presence of kinking all along the bloodline before use. In case the operator recognizes the presence of kinking, the tubings must be scrapped.

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**Actions Planned by Vital
Healthcare**

Re-train the workers on tubing rolling process as well as workers on 100% final inspection of finished products

**Action corrective de
Vital Healthcare**

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Some events of loose component inside primary packaging such as detached drainage bag and detached Transducer Protector (TP)

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Cause: Improper operation of the worker. Components may not be tightly screwed during production and become loose after long time transportation.

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Potential Risk: Contamination / Infection / Clotting

Risques potentiels : Contamination / infection/ coagulation

**Action to be Taken by
Customers / Users**

Operator must check the presence of detached components before open the primary packaging and take care during packaging opening that no loose components fall on contaminated surfaces (i.e floor).

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Actions Planned by Vital Healthcare

Retrain the workers. Increase the sampling amount for inspection.

Action corrective de Vital Healthcare

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Rare events of air entrance in the extracorporeal blood circuit. Majority parts of the events involve the transducer protector. The first TP was glued to the tubing and the second TP was screwed to the first TP manually. The second TP might not be firmly connected.

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Cause: Improper operation of the worker. Components may not be tightly screwed during production and become loose after transportation. If the connection is not tightened by operator according to section 5.4 of the manual, air entrance may occur.

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
Operator must check the tightening of all the bloodlines connection before use according to section 5.4 of relevant product *Instructions for Use*. The critical connection to check is illustrated below:

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**Actions Planned by Vital
Healthcare**

**Action corrective de Vital
Healthcare**

The double TP should be glued to ensure they're leak-proof and reduce the risk

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Rare events of not proper fitting between arterial chamber and machine holders.
Mauvaise tenue de la chambre artérielle dans son support.

Cause: There are different machine holder types for arterial chamber. One machine holder has several holes/places for fitting different chamber types. The arterial chamber may be placed at wrong hole/place of the machine holder and is not properly fitted in.

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Potential Risk: The falling-down of the chamber from the machine holder. The air

at the top of the chamber may come into the bloodline circuit.

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Action to be Taken by Customers / Users

Operator must check that the arterial chamber is placed at the right hole of the holder suiting its size and is properly fitted in.

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Actions Planned by Vital Healthcare

Develop an arterial chamber with same size as the ones that the users used previously.

Action corrective de Vital Healthcare

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Cause: The involved materials are harder than the lines which they used previously.

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Potential Risk: Difficulty to close the machine holder, malfunction on the air bubble detection.

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Action to be Taken by Customers / Users

Operator must check that the venous chamber is clamped/ fastened surely

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Causes:
Causes

- a) The double transducer protector is loosely screwed, causing air entrance in the tubing and clotting.

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circuit.

- b) If the blood volume level inside the arterial chamber and venous chamber is too high, the blood on top may get stagnant, which leads to clotting.

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- c) The clamp does not clamp tight, causing blood level inside the chamber to increase and the blood at surface may become stagnant, which may cause clotting.

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- d) NovaLine tubing's venous chamber filter is short filter, which may not allow too much blood pass through the filter at same time. When the blood volume is big at the venous chamber, the blood surface in the venous chamber may become clotted.

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- e) The drip chamber does not attach to the machine clamp well, causing the air alarm and clotting after pump stops.

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Adjust the anti-coagulation prescription if necessary. Monitor the blood volume inside the chamber, it is recommended the blood volume inside arterial chamber should be no more than 2/3 of the chamber volume and the blood volume/surface inside the venous chamber should be at same level as the dropper end. Once when coagulation is observed (by visual check or by an increasing venous pressure detected by monitor), stop the treatment, rinse back blood to patient in a safe manner and change the circuit according to indications provided by relevant *Instructions for Use*.

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Keep this Addendum with your NovaLine® Instructions for Use.

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stop the treatment rinse back blood to patient in a safe manner and change the circuit according to indications provided by relevant *Instructions for Use*.

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Please complete and return this form by FAX to the number below: +33 (0) 1 34 61 55 25 or by Email:
mv_france@baxter.com

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Subject: NovaLine® Tubing Sets for Hemodialysis – Instructions for handling problems encountered in clinical use

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