23 mars 2021

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The screenshot displays the 'User Defined Protocol Overview' table in the Ortho Clinical Diagnostics software. The table lists various tests with columns for Test, Original Test, Results, Last Modified, and Last Modifier. The row for 'UDP_4 ABO' is highlighted with a red circle. Below the table, the 'Modify' button in the bottom navigation bar is also circled in red.

Test	Original Test	Results	Last Modified	Last Modifier
UDP_Fya antitoxin	4 RAS Fya IgG	Fya	2/23/2021 11:32:20 AM	Admin SOL
UDP_Cw	08 RAS Lea Rvs	Anti-Cw	2/23/2021 4:34:46 PM	Admin SOL
UDP_08 RAS Fya IgG	08 RAS Fya IgG	Fya	2/23/2021 11:51:27 AM	Admin SOL
UDP_4 ABO(FWD/RVS)/R	4 ABO(FWD/RVS)/Rh-00	ABO Rh	2/23/2021 11:54:26 AM	Admin SOL
UDP_DS BV IAT IgG 08	DS BV IAT IgG 08 37 RRBC	DiSeries	2/23/2021 4:24:19 PM	Admin SOL
UDP_08 RAS Lea Rvs	08 RAS Lea Rvs	Lea	3/1/2021 5:00:35 PM	Admin SOL
UDP_4 ABO	4 ABO(FWD/RVS)/Rh-00	ABO Rh	3/3/2021 3:30:00 PM	Admin SOL

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The screenshot shows the 'Ortho Clinical Diagnostics' Admin SOL interface. The top navigation bar includes Home, Resources, Samples, Results, Errors, QC, Setup, Software, Maintenance, and Diagnostics. The user is logged in as 'Admin SOL' on '03/03/2021 15:32'. The main content area is divided into a left sidebar and a central workspace. The sidebar contains sections for 'Test Name*' (UDP_4 ABO), 'Assays*', 'Cassettes', and 'Result Interpretation'. The 'Assays*' section is highlighted with a red circle. The central workspace shows a 'Test Name' input field with 'UDP_4 ABO' and a virtual keyboard. The right sidebar contains a menu with 'General', 'Users', 'Testing' (highlighted), 'Results', 'System', 'Interfaces', 'Maintenance', and 'User Defined Reports'. The bottom navigation bar includes 'Back', 'Modify', 'Create Test', 'Delete Selected', 'Cancel', 'Save', 'Search', 'Help', and 'Resume Processing'.

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The screenshot shows the 'Ortho Clinical Diagnostics' Admin SOL interface. The top navigation bar is the same as in the previous screenshot. The user is logged in as 'Admin SOL' on '03/03/2021 15:33'. The main content area is divided into a left sidebar and a central workspace. The sidebar contains sections for 'Test Name*' (UDP_4 ABO), 'Assays*', 'Cassettes', and 'Result Interpretation'. The 'Assays*' section is expanded, showing a list of assay columns: 'Column 1*' (Anti-A), 'Column 2*' (Anti-B), 'Column 3*' (Anti-D), 'Column 4*' (Ctrl), 'Column 5*' (A1-Cells), and 'Column 6*' (B-Cells). The 'Column 5*' row is highlighted with a red circle. The central workspace is empty. The right sidebar contains the same menu as in the previous screenshot. The bottom navigation bar is the same as in the previous screenshot.

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The screenshot shows the Admin SOL interface for Ortho Clinical Diagnostics. The top navigation bar includes Home, Resources, Samples, Results, Errors, QC, Setup, Software, Maintenance, and Diagnostics. The 'Setup' menu is expanded, showing options like Profiles, Test Settings, Reagent Kits, QC Settings, and User Defined Protocols. The 'Test Settings' option is selected. On the right side, there is a sidebar with menu items: General, Users, Testing, Results, System, Interfaces, Maintenance, and User Defined Reports. At the bottom right, there are three buttons: Search, Help, and Stop Processing. The 'Stop Processing' button is circled in red.

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The screenshot shows the Admin SOL interface for Ortho Clinical Diagnostics. The top navigation bar is the same as in the previous screenshot. The 'Test Settings' option is selected in the Setup menu. The 'Test Settings' option is circled in red. On the right side, the sidebar menu is the same. At the bottom right, there are three buttons: Search, Help, and Resume Processing. The 'Show Details' button is circled in red.

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The screenshot displays the 'Ortho Clinical Diagnostics Admin SOL' interface. The top navigation bar includes icons for Home, Resources, Samples, Results, Errors, QC, Setup, Software, Maintenance, and Diagnostics. The main content area is divided into two panels: 'Automatic Result Acceptance' and 'Manual Result Acceptance'. The 'Automatic Result Acceptance' panel contains a table of settings, with the 'Require Manual Review of UDP' row circled in red. The 'Manual Result Acceptance' panel shows a single setting: 'Enforce Different Accept User' set to 'Yes'. A right-hand sidebar lists various system settings like General, Users, Testing, Results, System, Interfaces, Maintenance, and User Defined Reports. At the bottom, a navigation bar includes a 'Back' button, a circled 'Edit Automatic Result' button, and buttons for Search, Help, and Resume Processing.

Automatic Result Acceptance	
Enforce Manual Review (all results and cassettes)	No
Enforce Manual Review (all results and cassettes with error grades)	No
Auto Acceptance/Manual Review Rules (ignored if Enforce Manual Review is enabled)	
All Maintenance Tasks Must Be Completed Successfully	No
All Column Reactions Must Be Without Errors	Yes
Valid Test Result Interpretation (includes nondiscrepant columns)	Yes
Above/Below Positive Reaction Threshold	No
Enforce Different Accept User	No
Require Manual Review of UDP	No

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