Confirmation of correspondence between PCR methods: Boiling method and Non boiling method

Dear our valued customer,

We, Seegene Inc., hereby inform that the new extraction-free method which excludes 98 °C heating step for 3 min is compatible with the current extraction-free method.

In addition, while performing PCR, film can be used as an alternative of caps for PCR plates.

As shown below table, the new method leads to save hands-on time caused by boiling and chilling steps, and moreover it enhances Ct vales comparing to current method.

Therefore, the new extraction-free method is considered to be a way of increasing the efficiency of COVID test, in the situation where the demand of testing for COVID-19 has increased dramatically.

Eugene Kim, Ph. D

Managing Director of Marketing Strategy Department

Seegene Inc.

[Appendix] Method and Result

Current method	New method
1)Aliquot 45ul of Nuclease-free water into PCR tubes 2)Add 15ul of each sample into the tube containing aliquot of the Nuclease-free water 3)Close the cap, quick vortex and briefly centrifuge 4)Incubate at 98℃ for 3 min and Chill at 4℃ for 5 min 5)Prepare the Reaction Mastermix: use 5ul Sample's nucleic acid	1)Aliquot 45ul of Nuclease-free water into PCR tubes 2)Add 15ul of each sample into the tube containing aliquot of the Nuclease-free water 3)Prepare the Reaction Mastermix: use 5ul Sample's nucleic acid

Center	Positive concordance	positive thermolysis / Negative NO thermolysis	Negative thermolysis / positive NO thermolysis	Negative concordance	Total
#1	3	0	4	1	8
#2	61	1	0	11	73
#3	53	4	6	103	166
#4	35	1	11	47	94
#5	32	1	4	35	79
#6	15	0	3	28	46
#7	52	3	19	65	139
#8	53	2	3	66	124
#9	95	0	5	0	100
Total	399*	12**	55***	356	829

* Ct NO thermolysis < between 2~5 cycles versus thermolysis

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