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Genay, **October 9, 2019**

URGENT : FIELD SAFET NOTICE

This notice needs to be passed on all those who need to be aware within your organization or to any organization where the potentially affected devices have been transferred (if appropriate)

DUAL MOBILITY LINER diameter 28 size 48

Our Ref.: EN 020/F/ **269412800**

Dear Customer,

Details on affected devices HQN DM848 lot 269412800

This is to inform you of a product recall involving **HQN DM848 lot 269412800** manufactured by Groupe lépine.

Groupe Lépine has identified a potential head rotation defect (stainless steel or Alumina) after impaction. The investigation is underway; the main hypothesis is random machining issue. The defect is not reproducible so all the units in the batch are not affected by the anomaly; however we would like to carry out additional dimensional checks.

We also confirm that this defect is 100% detectable by the surgeon, the leaflet produced as the surgical technique requiring verification of the mobility of the head in the polyethylene core.

Furthermore, once mobility is obtained, there is no risk of blocking the head afterwards.

Finally, in the event of an anomaly is detected, the patient impact is very limited insofar as we systematically deliver for each surgery the same references in non-impacted or pre-impacted version, the surgeon being able to choose one or the other version in the event of an identified mobility problem.

As a precautionary measure and in order to carry out additional checks, Groupe Lépine has decided to recall lot 269412800.

Therefore, we kindly ask you:

- Identify and isolate the HQN DM848 devices lot 269412800 immediately upon receipt of this letter,
- Complete the product return form in the appendix and return it to us as soon as possible (by fax or mail). We will contact you to organize at our expense the return of these products and their replacement,
- Forward this letter to any other organization concerned by this recall (distributor, healthcare establishment, etc.).

Regarding patients who were implanted with a prosthesis from the recalled lot, Groupe Lépine in view of the above does not recommend any specific follow-up of the patients.

Pursuant to Article L1111-2 of the Public Health Code, Groupe Lépine invites the health professionals concerned to take this information into account and take the appropriate measures.

We also remind you of the need to report any undesirable effect observed with these devices to the National Agency for the Safety of Medicines and Health Products - Surveillance Department - by email to materiovigilance@ansm.sante.fr or by fax to 33 (0) 1 55 87 37 02.

Please accept our apologies for the inconvenience this may cause you and thank you for your confidence. Please accept, Madam, Sir, the expression of our highest consideration.

Laurence Fiscus

Quality and Regulatory affairs Director

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HQN DM848 lot 269412800

Please check ALL appropriate boxes.

- I have read and understand the recall instruction provided in the 'reference' letter
- I have checked my stock and have quarantined inventory consisting of:
 - No more stock
 - Quantity left:

| Reference | Lot number | Quantity |
|-----------|------------|----------|
| | | |
| | | |
| | | |
| | | |

- Indicate disposition of recalled product:
 - Returned
 - Destroyed
- I have identified and notified my customers that were shipped or may have been shipped this product.

Name: _____

Position: _____

Firm Name: _____

Address: _____

Please fax or mail completed return response form to:

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Laurence Fiscus – Vigilance officer

175 RUE JACQUARD – CS 50307 – 69727 GENAY CEDEX – FRANCE

Phone: +33 (0)4 72 33 02 95 – Fax: + 33 (0)4 72 35 96 50