04 juin 2020

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June 4, 2020

**URGENT: MEDICAL DEVICE FIELD SAFETY NOTICE**

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Dear Charlotte Charuet,

**Purpose of this Letter:**

TSI records indicate the products listed in Table 1 were distributed to you. The purpose of this letter is to advise you that TeDan Surgical Innovations, LLC (TSI) is voluntarily issuing a field safety notice for product exhibiting a potentially open pouch seal specific to the part and lot numbers indicated in Table 1.

*Table 1. Distributor Listing for Affected Part and Lot Number*

Item #	Product Name	Lot #	UDI #	Box Qty	Expiration Date
DS-0014	PHANTOM CS DISTRACTION SCREW, 14 MM, STERILE, 5/BOX	2018020801	00851797006450	30	1/31/2023
DS-0014	PHANTOM CS DISTRACTION SCREW, 14 MM, STERILE, 5/BOX	2018061101	00851797006450	30	5/31/2023

**Intended Use:**

TSI Distraction Screws are classified as Class I exempt devices according to 21CFR 888.4540; Orthopedic Manual Surgical Instrument (code: LXH). The TSI Distraction Screws are medical devices intended to be used for distraction in the anterior approach to the cervical spine. The screws are designed for use with TSI distractors and screwdrivers. These products are intended to be used by trained surgeons. The screw is exposed to the surgical area only and is removed prior to the end of the procedure. The screw is not left in the body.

TSI has had no reports of deaths or serious injuries or injuries related to this device.

**Reasons for the Field Safety Notification:**

The items identified for the field safety notification are labeled sterile. However, TSI has found that there is the potential for product pouches not to be sealed, potentially exposing the product to non-sterile conditions. This situation could result in a risk of infection to patients.

TSI is aware of 3 pouches out of 10,165 (potential product impacted resulting a 0.0295% probability) where there was no seal on the screw pouch opposite the chevron manufacturing pouch seal as illustrated in Figure 1. This issue was detected prior to use by the customer during visual package integrity inspection as required in the TSI Instructions for Use. TSI has not detected any lack of seals in any additional lots; however, the root cause was linked back to oversight by the contract packaging and

sterilization vendor for a specific time period with specific lots included in this field safety notice. As stated above, there have been no reports of deaths or serious injuries or injuries related to this device.

**Actions to be Taken by the Distributor:**

Specific to product in your inventory listing in Table 1, TSI requests that an immediate visual inspection of the pouch (Letter A shown in the red outline opposite the chevron seal [Letter B]) and look for an open seal. Any product exhibiting an open seal is to be quarantined, and then returned to TSI. Please contact TSI for RMA and shipping instructions.

It should be understood not all TSI Distraction Screws subject to this notice may have screw pouches with an open seal.

Devices with pouches that pass a visual inspection are sterile and may be used in accordance with the device instructions for use.

Specifically, TSI requests the following actions be performed to visually inspect the product packaging:

1. Carefully open the Distraction Screw box only on one (1) end, ensuring box integrity.
2. Remove the five (5) sterile pouches.
3. Visually inspect the seal zone (indicated as letter A in Figure 1) opposite the chevron seal (Letter B) for each sterile pouch.
4. IF the pouch is not sealed, **DO NOT USE THE DEVICE. Quarantine the device and contact TSI to receive an RMA, shipping instructions and replacements.**
5. IF the pouch is sealed, is not opened, the product remains sterile and may be used.

**Note: TSI will provide, by request following receipt of the acknowledgement, an additional "inspected by" label to reseal for those boxes where the all the pouches are sealed.**

6. Complete and return to TSI the enclosed FIELD SAFETY NOTICE ACKNOWLEDGMENT via fax or e-mail.

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
Figure 1. Drawing of Pouch Area to Inspect – Letter A

**Contact Information for Returns:**

<b>Name</b>	Account Management
<b>Telephone</b>	713-726-0886
<b>Email Address</b>	acctmgmt@tedansurgical.com
<b>Fax</b>	713-726-0846

In addition, this field safety notice should be carried out to the end user level. If you have distributed the identified product in Table 1 to end users, immediately identify all customers and notify them of this product field safety announcement. Please use the **attached customer letter**. Your assistance is appreciated and necessary to prevent any possible injury to patients. Also, please provide to these customers the enclosed acknowledgement response form for their completion and return to TSI.

Other TSI distraction screws not identified as part of this field safety notice may be used. However, as a precaution, always follow the instructions for use and carefully inspect product sterile packaging for any breaches as instructed in the device instructions for use. If breaches are identified DO NOT USE THE PRODUCT and contact TSI immediately.

**Other Information**

Please complete and return the enclosed mandatory acknowledgement and receipt form as soon as possible but not later than **July 6, 2020**.

If you have any questions, contact Lynne Davies, [ldavies@tedansurgical.com](mailto:ldavies@tedansurgical.com) or call 713-554-2130.

This field safety notice is being made with the knowledge of the Food and Drug Administration.

Name: (Print) Lynne Davies

Signature: *Lynne Davies*

Title: Vice President, Regulatory Affairs and Quality Assurance

**Enclosures:** TSI Acknowledgment and Receipt Form - 3005726885  
Customer Letter - 3005726885

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