| | |
|-----------------|---|
| EN 2 | Cook® EchoTip® Double Lumen Aspiration Needle Guidance for use, warnings and precautions |
| BG 3 | Cook® EchoTip® Double Lumen Aspiration Needle Указания за употреба, предупреждения и предпазни мерки |
| CS 5 | Cook® EchoTip® Double Lumen Aspiration Needle Návod k použití, varování a bezpečnostní opatření |
| DA 6 | Cook® EchoTip® Double Lumen Aspiration Needle Brugsvejledning, advarsler og forholdsregler |
| DE 7 | Cook® EchoTip® Double Lumen Aspiration Needle Anwendungsanleitung, Warnhinweise und Vorsichtsmaßnahmen |
| EL 9 | Cook® EchoTip® Double Lumen Aspiration Needle Οδηγίες χρήσης, προειδοποιήσεις και προφυλάξεις |
| ES 11 | Cook® EchoTip® Double Lumen Aspiration Needle Guía de uso, advertencias y precauciones |
| ET 12 | Cook® EchoTip® Double Lumen Aspiration Needle Kasutusjuhised, hoiatused ja ettevaatusabinõud |
| FI 13 | Cook® EchoTip® Double Lumen Aspiration Needle Käyttöohje, varoitukset ja varotoimet |
| FR 14 | Cook® EchoTip® Double Lumen Aspiration Needle Directives d'utilisation, avertissements et mises en garde |
| HR 15 | Cook® EchoTip® Double Lumen Aspiration Needle Priručnik za uporabu, upozorenja i mjere opreza |
| HU 16 | Cook® EchoTip® Double Lumen Aspiration Needle Útmutató a használatához, a figyelemztetésekhez és az óvintézkedésekhez |
| IT 17 | Cook® EchoTip® Double Lumen Aspiration Needle Linee guida per l'uso, avvertenze e precauzioni |
| LT 18 | Cook® EchoTip® Double Lumen Aspiration Needle Naudojimo nurodymai, įspėjimai ir atsargumo priemonės |
| LV 19 | Cook® EchoTip® Double Lumen Aspiration Needle Lietošanas norādījumi, brīdinājumi un piesardzības pasākumi |
| NL 20 | Cook® EchoTip® Double Lumen Aspiration Needle Aanwijzingen voor het gebruik, waarschuwingen en voorzorgsmaatregelen |
| NO 21 | Cook® EchoTip® Double Lumen Aspiration Needle Veiledning for bruk, advarsler og forholdsregler |
| PL 22 | Cook® EchoTip® Double Lumen Aspiration Needle Wskazówki dotyczące użycia, ostrzeżenia i środki ostrożności |
| PT 23 | Cook® EchoTip® Double Lumen Aspiration Needle Orientações de utilização, advertências e precauções |
| RO 24 | Cook® EchoTip® Double Lumen Aspiration Needle Instrucțiuni de utilizare, avertismente și precauții |
| RU 25 | Cook® EchoTip® Double Lumen Aspiration Needle Руководство по эксплуатации, предупреждения и меры предосторожности |
| SK 27 | Cook® EchoTip® Double Lumen Aspiration Needle Pokyny na použitie, varovania a bezpečnostné opatrenia |
| SL 28 | Cook® EchoTip® Double Lumen Aspiration Needle Smernice za uporabo, opozorila in previdnostni ukrepi |
| SV 29 | Cook® EchoTip® Double Lumen Aspiration Needle Bruksanvisning, varningar och försiktighetsåtgärder |
| TR 30 | Cook® EchoTip® Double Lumen Aspiration Needle Kullanım, uyarılar ve önlemler için kılavuz |

COOK® ECHOTIP® DOUBLE LUMEN ASPIRATION NEEDLE

1 DEVICE DESCRIPTION

The Cook® EchoTip® Double Lumen Aspiration Needle consists of: a stainless-steel double lumen needle with a bevelled tip for puncture and an echogenic tip for enhanced visualisation under ultrasound, a manipulation handle, aspiration line, a silicone stopper with Luer fitting to allow connection to a vacuum source and flushing line.

2 INTENDED USE

The Double Lumen Ovum Pick-up Needles are intended to be used for ultrasound guided transvaginal simultaneous or independent aspiration and flushing of ovarian follicles during ovum pick-up procedures.

3 CONTRAINDICATIONS

The Double Lumen Ovum Pick-Up Needles are not for use on patients:

- With an active vaginal or intrauterine infection
- With a sexually transmitted disease
- With a recent uterine perforation
- With a recent cesarean section
- After recent pelvic inflammatory disease
- Who are currently pregnant

4 WARNINGS

- The Double Lumen Ovum Pick-Up Needles should only be used by persons or under the guidance of qualified persons having adequate training and familiarity with ovum pick-up procedures.
- There is not sufficient clinical evidence of safety and/or efficacy to support the use of Double Lumen Ovum Pick-Up Needles for transabdominal ovum aspiration procedures.
- There is not sufficient clinical evidence of safety and/or efficacy to support the use of Double Lumen Ovum Pick-Up Needles for independent cystic aspiration procedures.
- Cases of significant hemoperitoneum post-operatively have rarely occurred hours to days after oocyte retrieval. All patients should be screened pre-operatively for increased risk of bleeding. At completion of the procedure a final ultrasound review of the ovaries and the recto-uterine pouch should ensure no active bleeding is seen. All patients with increased pain post-operatively should be reviewed.
- If an endometrioma or a hemorrhagic follicle is inadvertently punctured, the needle and collective tubing should be immediately withdrawn and flushed, and the collection tube should be changed. Treat patients in accordance with the higher risk of infection.
- The enlarged ovary following the controlled hyperstimulation as part of infertility treatments is at risk of twisting and restricting the ovarian blood supply. Tangential needle insertion may increase the risk of that ovarian torsion. Firm pressure on the ovary with the ultrasound transducer may assist to stabilise the ovary. To further stabilise a mobile ovary in one place supra-pubic pressure can be applied (with the help of an assistant).
- This single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to device failure and/or infection.

5 PRECAUTIONS

- Hematuria may occur due to the needle penetrating a filled bladder during transvaginal ultrasound aspiration. This complication typically resolves spontaneously within a day.
- Extravasation of urine may occur within the abdominal cavity if a needle puncture traverses the bladder. Patients should be monitored for evidence of this known complication; however, there is typically no prolonged associated discomfort or adverse sequelae.
- Infection may be introduced via needle puncture and result in urinary tract infection (UTI), pelvic inflammatory disease (PID), peritonitis, uterine infection or cystitis. Recommendations to minimise occurrence include use of only biocompatible materials, flushing the needle prior to use (and any other accessories used) with compatible sterile culture media, and closely adhering to sterile techniques.
- Vaginal/cervical bleeding has been reported to be associated with the transvaginal route for oocyte retrieval via needle aspiration. Bleeding is typically easily controlled with direct pressure or by a suture.
- Pelvic hematomata have been reported in the vaginal vault within an ovary, under the bladder and in the broad ligament. Patients should be monitored for evidence of this known complication; however, there is typically no further intervention required.
- Due to the larger gauge of double lumen cannula, additional intra-operative pain management may be required.
- Ultrasound equipment must be of sufficient resolution to visualise the needle tip and pelvic anatomy. If the needle tip and pelvic anatomy cannot be visualised the procedure cannot be successful.

6 HOW SUPPLIED

STERILE – DO NOT RESTERILIZE – SINGLE USE ONLY

Double Lumen Ovum Pick-Up Needles are supplied sterile in peel-open packages. They have been sterilised using Ethylene oxide and are intended for single use only.

Do not resterilise the device. Do not use after the expiration date printed on the label. Store in a dark, dry, cool place. Avoid extended exposure to light and extremes of temperature and humidity.

7 INSPECTION OF DEVICE

1. Check that the device packaging is not damaged and that the 'use by' date has not expired.
 - CAUTION: Do not use if the sterile packaging is damaged or unintentionally opened before use.**
 - CAUTION: Do not use after the expiration date printed on the label.**
2. Carefully remove the needle from the packaging maintaining the sterility of the product.
3. The sterile needle should be inspected for irregularities in needle tip profile or kinking of any supplied tubing.

8 DEVICE PREPARATION

1. Connect the Luer fitting on the silicone stopper to a new disposable vacuum line with hydrophobic filter attached to vacuum pump.
2. Ensure the vacuum pump has been prepared as per respective manufacturer's instructions.
3. Fit the silicone stopper onto the collection tube (designed to fit 14 mL Falcon tubes).
4. Place the needle tip in a spare collection tube containing approximately 5 mL of sterile culture medium and apply vacuum. Ensure culture medium is aspirated into collection tube.

CAUTION: Ensure all connections are tight and leak proof. Failure to confirm connections and system patency may result in failure to complete the procedure.

NOTE: The vacuum pump pressure should be calibrated to give the desired flow rate for the specific type of needle being used as per the vacuum pump manufacturer's instructions.

NOTE: If blockages occur, remove silicone stopper from collection tube and connect a luer syringe filled with sterile culture medium into the aspiration line port. Inject the culture medium and attempt to force the blockage out. Refer to the vacuum pump manufacturer's instructions for additional blockage troubleshooting. If blockage persists, remove blocked device and replace with new device for procedure completion.

5. Discard the used collection tube and replace silicone stopper into collection tube primed with sterile culture medium.
6. Ensure that the ultrasound equipment, ultrasound gel, transducer protective sheath and needle guide have been prepared as per respective manufacturer's instructions.

9 DEVICE OPERATION

1. Position the patient in the lithotomy position on operating table. Local or general anaesthetic may be administered as necessary.
2. Introduce an ultrasound transducer into the vaginal fornix to visualize the ovary and follicles. Identify the follicles to be aspirated. Check for the presence of blood vessels in and around the ovary and determine a direct path into the ovarian follicles to be aspirated.
 - NOTE:** For mobile ovaries, firm pressure with the ultrasound transducer may assist to stabilize the ovary.
3. Check the needle guide is appropriately placed on the ultrasound transducer with no gaps and that the transducer protective sheath does not hinder path of the needle projection.
 - CAUTION: In passage of the needle between personnel, be aware of needle tip location to avoid needle stick injury.**
4. Insert the needle into the needle guide, ensuring there is smooth movement along its length.
5. Visualize the ovary and follicles to be aspirated and advance the needle tip to puncture the vaginal wall.
 - WARNING: The needle tip must be visible at all times as it is maneuvered to the ovaries and into each follicle. If visibility is lost, do not advance the device further. Manipulate the ultrasound transducer to obtain different views. If still unable to visualize the needle tip, withdraw the needle from the patient.**
 - NOTE:** The combination of the needle bevel and EchoTip® enhances visualization of the position of the needle tip. The handle indent indicates bevel orientation as well as facilitating grip.
6. Line up the target follicle using the needle guide on the ultrasound monitor and advance the needle tip into the centre of an ovarian follicle via a rapid, stabbing motion.
 - WARNING: To further stabilize a mobile ovary in one place, supra -pubic pressure can be applied (with the help of an assistant).**
 - WARNING: Be aware that big iliac vessels may resemble a follicle on an ultrasound image when in a cross-sectional plane. Rotate the ultrasound transducer to visualize anatomy to avoid inadvertent puncture.**
7. Apply vacuum to aspirate the follicular contents into the collection tube. As the follicle collapses, rotate the needle tip within the follicle to ensure complete emptying of the follicle.
 - NOTE:** To avoid aspiration difficulty, ensure tubing does not become kinked.
 - NOTE:** The vacuum pump pressure should be calibrated to give the desired flow rate for the specific type of needle being used as per the vacuum pump manufacturer's instructions.
 - NOTE:** If blockages occur, stop aspiration, remove the needle from the patient, remove silicone stopper from collection tube and connect a luer syringe filled with sterile culture medium into the aspiration line port. Inject the culture medium and attempt to force the blockage out. This ejected fluid should be checked for oocytes. Refer to the vacuum pump manufacturer's instructions for needle blockage troubleshooting. If blockage persists, remove blocked device and replace with new device for procedure completion.
8. Move needle tip to the adjacent follicle and repeat **Step 7**, maintaining vacuum aspiration. Where possible, keep the needle tip within the stroma or follicles of that ovary until follicular aspiration is complete.
 - CAUTION: Avoid recurrent penetration of the vaginal vault and ovarian capsule to minimise risks of infection and bleeding.**
9. Replace the silicone stopper into primed collection tubes at appropriate aspiration intervals. To avoid aspiration difficulty, ensure aspirated contents do not reach the lip of silicone stopper in collection tube.
10. If required the follicle can be flushed, as described below.
 - a. Stop aspiration.
 - b. Use a Luer syringe filled with a sterile follicle-flushing buffer. Insert the syringe into the Luer fitting on the flushing line.
 - c. With the needle tip in the collapsed follicle, slowly inject (1-2 mL per second) the flushing buffer to refill the follicle.
 - d. Aspirate the follicular contents.
11. Confirm all follicles in the ovary have been aspirated. Stop aspiration before withdrawing the needle through the vaginal vault.
12. Place the needle tip in a spare collection tube containing approximately 5 mL of sterile culture medium and apply vacuum. Ensure culture medium is aspirated into the collection tube.
 - NOTE:** This aspirated content should be checked for oocytes.
13. Reposition the ultrasound transducer to visualize the remaining ovary. Repeat **Steps 5 to 10**.
14. At the completion of the aspiration procedure, stop aspiration and remove the needle from the patient.
15. Place the needle tip in a spare collection tube containing approximately 5 mL of sterile culture medium and apply vacuum. Ensure culture medium is aspirated into collection tube.
 - NOTE:** This aspirated content should be checked for oocytes.
16. Before the ultrasound transducer is removed, check for any evidence of continuing bleeding in the ovaries or peritoneal cavity.
17. Check the vagina for signs of bleeding. If vaginal bleeding is experienced post procedure, apply steady pressure. Consider inserting a suture to control continuing vaginal vault or cervical bleeding.

10 DISPOSAL OF DEVICES

1. Dispose of aspiration needle in appropriate sharps container and any consumables in accordance with institutional guidelines.
 - WARNING:** This single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to device failure and /or infection.

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