	Field Safety Notice	FSN Reference number
		22-001 rev. A
		Page 1 sur 4

AB MEDICA

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Our records indicate that your facility has received one or more products affected by this safety advisory notice. The table below provides a complete list of affected products.

Produit concerné / Concerned product	Référence / Reference	Numéro de lot / batch number
Insert Crochet pneumo-dissection D=3,5mm Lg=33cm plus tube, sterile Usage unique Hook insert and tube for pneumodissection D=3,5mm Lg=33cm steril single use	C35-ID-CH-33-ST	EK07210004
Insert Crochet pneumo-dissection D=5mm Lg=33cm plus tube, sterile Usage unique Hook insert and tube for pneumodissection D=5mm Lg=33cm steril single use	C05-ID-CH-33-ST	EK07210006

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During the qualification of the aging, AB Medica sas identified a non-conformity related to a European standard (NF EN 868-5 annex E(2018) of the sealing width of its blister/Tyvek packaging

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An initial investigation was conducted and concluded that the TYVEK batch used for this validation did not meet AB Medica's specifications. A second investigation determined that production batches were also affected by this non-conformity because the TYVEK used belonged to the same batch.

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Internal investigations have shown by toluidin blue test that the packaging is still sealed despite this non-compliant sealing width. More, resistance of sealing tests have been performed proving the strength of the sealing.

	Field Safety Notice	FSN Reference number
		22-001 rev. A
		Page 2 sur 4

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The product instructions for use (NU-PNE-EN_G) contains the following statement: "Check the integrity of the sterile packaging before use. Do not use if sterile packaging is damaged". The symbol  appears on the product packaging: Do not use if package is damaged.

	Field Safety Notice	FSN Reference number
		22-001 rev. A
		Page 3 sur 4

Actions requises / Required actions :

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In order to minimise the impact on patients we ask that the following measures be taken:

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This notice needs to be passed on all those who need to be aware within your organization or to any organization where the potentially affected devices have been transferred.

Please transfer this notice to other organizations on which this action has an impact.

Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action. (if appropriate)

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Thank you in advance for your prompt attention to this matter. Please contact your company correspondent within 5 business days to confirm quantity at your location. If you have questions related to this issue, please contact Vincent LERE at lere.vincent@abmedica.fr or +33 (0) 2 48 52 59 57

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Please return the products to the following address:

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5. Le cas échéant, des dispositifs de remplacement exempts de ce problème vous seront fournis par AB medica.
If necessary, AB medica will provide you replacement devices free of this problem.
AB MEDICA SAS