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new case, keep base data

Version 2.7en
2012-12-03

1 Administrative information
To which NCA(s) is this report being sent? ANSM(FR), BSI(nb), MHRA(GB), AFMPS(BE), BFARM(DE), AEMPS(ES), SANITA(IT), FIMEA(FI), SMPA(SE), SwissMedic (CH)
Type of report ☒ Initial report ☐ Follow-up report ☐ Final report
Date of this report 2022-07-06
Reference number assigned by the manufacturer VR22-1
FSCA reference number assigned by NCA
Incidence reference number assigned by NCA
Name of the co-ordinating NCA Competent Authority (if applicable)

2 Information on submitter of the report
Status of submitter ☒ Manufacturer ☐ Authorised Representative within EEA and Switzerland ☐ Others: (identify the role)

3 Manufacturer information		new
Name Renishaw Mayfield SARL		
Contact Name Stéphane Vinot		
Address 31, rue Ampère		
Postcode 69680	City Chassieu	
Phone +33(0)4 37238900	Fax	
E-mail stephane.vinot@renishawmayfield.com	Country FR - France	

4 Authorised Representative Information

new

Name	
Contact Name	
Address	
Postcode	City
Phone	Fax
E-mail	Country DE - Germany

5 National contact point information

new

National contact point name Renishaw Mayfield SARL	
Name of the contact person Stéphane Vinot	
Address 31, rue Ampère	
Postcode 69680	City Chassieu
Phone +33(0)4 37238900	Fax
E-mail stephane.vinot@renishawmayfield.com	Country FR - France

6 Medical device information

new

Class

- ☐ AIMD Active implants
- ☐ MDD Class III
- ☒ MDD Class IIb
- ☐ MDD Class IIa
- ☐ MDD Class I
- ☐ IVD Annex II List A
- ☐ IVD Annex II List B
- ☐ IVD Devices for self-testing
- ☐ IVD General

Nomenclature system (preferable GMDN)

GMDN

Nomenclature code

32568

Nomenclature text

Neurological Stereotactic Surgery System

Commercial name/ brand name / make

neuroinspire (accessory of the neuromate system)

Model number

N/A

Catalogue number

N/A

Serial number(s)

N/A

Lot/batch number(s)

N/A

Device Mfr Date**Expiry date****Notified Body (NB) ID-number**

BSI 2797 (for Europe) & BSI 0086 (for GB)

Accessories / associated devices (if applicable)

neuroinspire can be supplied as non robot system. It is also incorporated into the stereotactic system neuro|mate robot and is placed onto the market as an integrated system by Renishaw Mayfield.

Software version number (if applicable)

All software versions

7 Description of the FSCA	
Background information and reason for the FSCA When an enhanced DICOM format sagittal image series is used within surgical planning, the series is flipped (right/left) resulting in the image becoming a mirror of the actual image. Within surgical planning and image fusion this will cause unacceptable fusion with the other image series used for the planning, due to the incorrectly displayed image series. The flipped image will not line up with the other image series and will not provide a clinically acceptable 3D image of the head/brain within the stereotactic space. The clinician is informed of the fusion and is required to accept the fusion result based on their clinical judgment. The workflow allows the images and fusion to be reviewed further and image series to be removed and/or exchanged as required. Unacceptable fusion with the other image data series would be expected due to the non-symmetrical nature of the human head and brain. It is not possible to create a plan with sagittal image series only as this would not allow positioning of the patient in relation to the frame/robot space. Unacceptable fusion would result in a delay to the surgery as the clinician obtains an additional scan series and completes a new fusion of the images. Only if the fusion, using an enhanced sagittal image series, is accepted and used in the surgery, could the implantable be delivered to a non-intended target, this may result in death or serious deterioration of health.	
Description and justification of the action (corrective / preventive) An FSN will be provided to all users to explain the situation and advise against using enhanced DICOM format series within the software. This action will prevent the potential hazard from occurring. A root cause investigation will be undertaken to address the issue within the software and identify a solution.	
Advice on actions to be taken by the distributor and the user The FSN provides the below actions that are to be taken by the customer. Images stored in enhanced DICOM format should not be imported into neuroinspire until a solution has been implemented. As per the current IFU, care should continue to be taken when reviewing the outcome of the fusion process between image sets to ensure that the solution is clinically suitable.	
Progress of FSCA , together with reconciliation data (Mandatory for a Final FSCA)	
Time schedule for the implementation of the different actions FSN to be sent to all Renishaw Mayfield customers by 31 July 2022 (Will be sent before this date depending on the time needed for the translation) A fix is expected to be installed before end of December 2022	
Attached please find ☒ Field Safety Notice (FSN) in English ☐ FSN in national language ☐ Others (please specify)	FSN Status ☐ Draft FSN ☒ Final FSN

The medical device has been distributed to the following countries:

within the EEA and Switzerland

☐ AT	☒ BE	☐ BG	☒ CH	☐ CY	☐ CZ	☒ DE	☐ DK
☐ EE	☒ ES	☒ FI	☒ FR	☒ GB	☐ GR	☐ HU	☐ IE
☐ IS	☒ IT	☐ LI	☐ LT	☐ LU	☐ LV	☐ MT	☐ NL
☐ NO	☐ PL	☐ PT	☐ RO	☒ SE	☐ SI	☐ SK	☐ TR

Candidate Countries

☐ HR

☐ All EEA, candidate countries and Switzerland

Others:

8 Comments

Submission of this report does not, in itself, represent a conclusion by the manufacturer and/or authorised representative or the National Competent Authority that the content of this report is complete or accurate, that the medical device(s) listed failed in any manner and/or that the medical device(s) caused or contributed to the alleged death or deterioration in the state of the health of any person.

Signature

I affirm that the information given above is correct
to the best of my knowledge

print

check

sent on Wednesday, July
06, 2022 18:39:46