****FIELD SAFETY NOTICE****

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Patrick QUERO
Quality & Regulatory Affairs Director
AL.CHI.MI.A. Srl

24 JUILLET 2024

****FIELD SAFETY NOTICE****

****Identification of the affected product:****

Product name: TISSUE-C
Product reference: TIS 001-00
FSN reference : FSN_07_2024_TIS 001-00

****RECIPIENT:****

****Date:****

24/07/2024

Dear Users,

We are sending you this notification to inform you of a heterogeneity of results observed when using the product TISSUE-C in FRANCE, referenced as [TIS 001-00].

**Problem description: **

It has been identified that some users may experience significant variability in results in terms of corneal transplantation rates when using the TISSUE-C product.

This heterogeneity appears to be related not to product quality issues, but to variations in product implementation procedures and variable quality controls in terms of test type, time and methods applied. Incorrect use of the product and non-standardized quality controls can lead to cell loss in tissues (e.g., when used as a transport medium for TISSUE-C), thus compromising the reliability and repeatability of the results obtained.

**Recommendations for correct use: **

To minimize the variability of results and ensure correct use of the product, we ask you to carefully follow the following recommendations:

1. ****Review the instructions for use: **** Make sure that all procedures for using the TISSUE-C product are followed in accordance with the instructions provided in the user manual ref.: MCM246 rev.10 (storage temperature **between -25°C and -15°C** and culture temperature 31°C).
2. ****Preparation of the product:**** Ensure that the TISSUE-C product is thawed at **room temperature** and/or in the refrigerator (**2-8°C**). Shake before use.
3. ****Medium staining: **** Do not use if the solution is yellow, purple or cloudy.
4. ****Use of the product: **** Ensure that the tissue is kept out of contact with the bottom of the vial during the culture stage. This is essential for the proper conservation and preservation of the tissue. The *CornealFloat* device, included with the TISSUE-C product, has been specially designed for this step.
5. ****Temperature : **** Respect the temperature indicated by the manufacturer during the Organoculture phase: **31°C**.
6. ****Optimise quality control: **** Rinse tissues with BSS and/or PSS-L Solution ref. GRS 003-00. In the event of reduced visualisation of endothelial cells during microscopic analysis, apply a hypotonic solution to improve visualisation of endothelial cells. Assess endothelial cell density and endothelial mortality during tissue storage using a low concentration Trypan Blue solution (**max 0.25%**) and a maximum application time of **30 sec** per analysis.
7. ****Sampling : **** Make sure that the volume of medium used during storage of the cornea (collection of samples for microbiological analysis) is greater than or equal to 80% of the initial volume (100ml). Transfer the cornea to a new TISSUE-C medium if the residual volume **is less than 80 ml**.

Page 3 of 5

8. **Transport:** Make sure to use CARRY-C medium to guarantee perfect preservation of the cornea before transplantation.
9. **Train users:** Organise training sessions for all potential users to ensure correct understanding and implementation of the above instructions.
10. **References :** Consult the following guidelines for further information.
 - a. <https://www.eeba.eu/minimum-medical-standards-revision-5>
 - b. <https://www.eeba.eu/technical-guidelines-for-ocular-tissue-revision-11>

Corrective actions taken:

We undertake to review the instructions for use in order to clarify good practice in the use of our products. We plan to provide some information in the notice next September.

Contacts for assistance:

If you have any questions or require further assistance, please contact our customer service department at +39 049 8962074 or by email at vigilance@alchimiasrl.com.

We thank you for your cooperation and apologize for any inconvenience caused. The safety and satisfaction of our customers is our top priority, and we are taking all necessary steps to resolve this problem as quickly as possible.

Users affected by this safety action have received the letter by post with acknowledgement of receipt to be returned to us duly completed.

NOTE: the ANSM has been informed of this safety action.

Thank you for your confidence in our products. Yours faithfully

Patrick QUERO
Quality & Regulatory Affairs Director
AL.CHI.MI.A. Srl

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