Date: 2025-06-19

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English version

Purpose:

This technical note explains how to clean and disinfect endoscope simulators for leak testing (SL-OT-L), and how to store them after each use.

I. Preparation:

Before use:

Swab or brush: For cleaning the inside of the leak test connector.

70% ethyl alcohol: To disinfect the leak-test connector (spray twice).

Use one of the following wipes to wipe the simulator's connectors and ensure bactericidal and yeasticidal activity:

- Incidin Oxy Wipe S
- Alcohol-based wipes (alcohol concentration above 50%)
- Wipes on oxidizing base (H2O2, chlorine)

After use:

Single-use dry wipe: to dry the external surface of the simulator.

Medical air: to remove moisture from inside the simulator tube.

To clean the simulator's socket after use, use one of the following wipes:

- Incidin Oxy Wipe S
- Alcohol-based wipes (alcohol concentration above 50%)
- Wipes on oxidizing base (H2O2, chlorine)
- Clean packaging: To keep the simulator dry.

II. Before use:

➤ Checking and cleaning the leak test socket:

Apply hygienic hand rub with a hydroalcoholic product. Check that the leak test socket (see figure1) is clean. If necessary, clean the inside of the connector with a swab or brush to remove any soil or debris.

Spray 70% ethyl alcohol (spray twice) directly onto the leak test connector from an approximative distance of 10 cm.

Be sure to cover all internal surfaces of the leak test connector with ethyl alcohol solution (see leak test socket disinfection protocol in the user manual).

➤ Contact time

Leave for 5 min to ensure effective action. Depending on the disinfectant used, please follow the contact time recommended by the manufacturer.

Document #: TEM-02567 Revision: 4

*Figure 1: Leak test socket***☒ Simulator tube cleaning:**

Before connecting the simulator tube, clean it using a recommended wipe (see paragraph I. Preparation) to wipe the simulator's connectors, making sure to cover the entire surface. Leave for the contact time required for bactericidal and yeasticidal activity.

Important reminder:

Make sure the inside of the simulator tube is dry before connecting it. In case of moisture, replace it with a new simulator.

➤ Contact time

Leave for 5 min to ensure effective action. Depending on the disinfectant used, please follow the contact time recommended by the manufacturer.

*Figure 2: Drop of water inside the simulator tube**Figure 3: Clean simulator tube with no water inside*

Document #: TEM-02567 Revision: 4

III. After use:**Removing and drying the simulator tube :**

Disinfect your hands hygienically with a hydroalcoholic product. Remove the simulator from the leak test connector in the bowl of the washer.

Dry the outer surface of the simulator with a dry disposable wipe and medical air. Make sure the inside of the simulator tube is dry. If moisture remains, replace it with a new simulator tube.

Cleaning the simulator connectors:

Use a recommended wipe (see paragraph I. Preparation) to clean the simulator's connectors, turning them around to apply a mechanical action and cover the entire surface with the disinfectant solution. Leave for the contact time required for bactericidal and yeasticidal activity.

Simulator tube storage:

Keep the simulator dry and in clean packaging.

Check that the leak test connector is clean. If necessary, clean the inside of the connector with a swab or brush to remove any soil or debris.

Spray 70% ethyl alcohol (spray twice) directly onto the leak test connector from an approximative distance of 10 cm.

Be sure to cover all internal surfaces of the leak test connector with ethyl alcohol solution (see leak test socket disinfection protocol in the user manual).

➤ WD SELF-DISINFECTION

Connect SL-OPE plug to leak test socket on bowl and launch cycle 6.

Writer	Reviewer	Approver
Name: M'Bark HAIMAD Date: 30/06/2025 Signature:	Name: Nacima NECHAK Date: 30/06/2025 Signature:	Name: Yolande MARTINEZ Date: 30/06/2025 Signature:

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ENGLISH VERSION

1. OBJECT

This technical note explains how to reprocess and keep the leak test socket in the bowl clean. It also describes the daily maintenance actions to be carried out. This will prevent any source of contamination and ensure optimal hygiene.

Reminder: Failure to comply with the maintenance instructions may invalidate the manufacturer's warranty and can adversely affect the quality of the device reprocessing.

2. DAILY MAINTENANCE

Daily maintenance is required to ensure optimal performance of Soluscope washer disinfectors. These actions must also be performed after each maintenance intervention on the machine.

- Wipe the outside of the washer disinfectant with a pH-neutral cleaning and disinfecting wipe.
- Clean the washer disinfectant plate, the interior edge of the bowl (first 5 centimeters from the edge) and the lid, using a surface disinfectant spray (e.g. Aniospray Surf 29 by Laboratoires Anios).
- Reprocess the leak test socket (*figure 1*) in the bowl (see reprocessing of the bowl leak test procedure paragraph 3).
- Clean the lid seal using a surface disinfectant spray.
- Check that there is no water in the leak test connector. Use a clean cloth soaked in 70% ethyl alcohol to disinfect the SL-OTE connector tip (*figure 2*). Make sure the connector is completely dry before plugging it in.
- Launch self-disinfection of the washer-disinfectant (cycle 6).
- Dry and check the SL-OTE connectors (*figure 2*), then disinfect them using a surface disinfectant spray. Make sure there are no cracks or kinking marks. Store the connectors in a clean and dry location. Lubricate the seal of leak test connectors with silicone oil (Lubrispray P recommended).

3. METHOD

3.1. PREPARATION

Wear gloves to protect your hands.

Make sure the 70% ethyl alcohol is in a ready-to-use spray bottle (not diluted in-house).

Disconnect all leak test connectors.

Important reminder: Using a water-based disinfectant as an alternative is strictly prohibited. It may damage the endoscope.

3.2. INITIAL CLEANING

If the leak test socket is visibly soiled, be sure to clean the inside of the socket with a swab or brush to remove any dirt or debris.

Figure 1: Leak test socket in the bowl

3.3. SPRAYING

Spray (2 sprays) 70% ethyl alcohol directly onto the leak test socket from approximately 10cm. Be sure to cover all internal surfaces of the leak test socket.

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3.5. BLOWING

a) For software versions up to v1.8.4

Do not connect anything on leak test socket in the bowl. Close the lid and launch cycle 11. After 2 minutes of blowing, an alarm 1 should appear. Return to main view.

b) For software versions v1.9.0 and above

Use cycle 10 to blow off any excess alcohol and wait for the end of cycle. At this step, wipe the lid seal.

3.6. WIPING

Use a clean cloth soaked in 70% ethyl alcohol to reprocess the SL-OPE plug tip.

Make sure the plug is completely dry before plugging it in.


Figure 2: SL-O TE area to be reprocessed

3.7. PRECAUTIONS

Avoid using too much alcohol to avoid flooding the leak test socket (Do not spray on contact with leak test socket).

3.8. WDE SELF-DISINFECTION

Connect SL-O TE plug to leak test socket on bowl and launch cycle 6.

Writer	Reviewer	Approver
Name: M. HAIMAD Date: 30/06/2025 Signature:	Name: N. NECHAK Date: 30/06/2025 Signature:	Name: Y. MARTINEZ Date: 30/06/2025 Signature: