FSN Ref. : FSN-2025-001 FSCA Ref. : FSN-2025-001

Date : 21 août 2025

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Date: 21 August 2025

Urgent Field Safety Notice (FSN)

To the attention of laboratory managers, local vigilance correspondents and directors of the health establishment

Device Type: Antimicrobial Susceptibility Discs

Commercial Name: Imipenem 10 µg IPM10

Device Model/Catalogue ref. number: 065048L

Manufacturer: Liofilchem

Distributor: INTELLIGENCE ARTIFICIELLE APPLICATIONS (i2a)

Dear Customer,

This letter is intended for all **Imipenem 10 µg IPM10 (ref. 065048L)** users. Our records indicate that your laboratory received the following products (reference and lot numbers included in the table below).

Table 1.

Ref.	Product Description	Lot No.	Expiry Date
065048L	Imipenem 10 µg IPM10 (SirScan Discs)	112023085	2026.11.19
		112023104	2026.11.19
		030124087	2027.03.01
		030424047	2027.03.04
		030424048	2027.03.04
		061224044	2027.06.12
		061224072	2027.06.12
		091924022	2027.09.19
		120424019	2027.12.04
		120424097	2027.12.04
		120624055	2027.12.06
		120624056	2027.12.06
		120924065	2027.12.09
		120924138	2027.12.09

Description of the issue

QC failures (zone diameter result below the QC range) for *Escherichia coli* ATCC 25922 and *Pseudomonas aeruginosa* ATCC 27853 on imipenem 10 µg IPM10 discs (ref. 065048L) reported from the field. A stability issue was observed for discs used within 7 days after the first opening (in-use shelf life).

Liofilchem conducted an internal investigation to confirm product issue and determine root cause.

The following have been identified:

- The investigation confirmed a potential performance issue for QC strains and clinical isolates on Imipenem 10 µg IPM10 (ref. 065048L, lots listed in Table 1).
- This issue is currently suspected to be caused by an inhomogeneity in the product, leading to faster rates of degradation of the antibiotic in some discs while others are not affected.

Clinical impact

There is a potential to obtain a zone diameter that is smaller than expected during the product's in-use shelf-life. This type of error would be detectable during quality control testing as an out of specification zones would be obtained. Patient results may also be affected leading to erroneous categorization (resistant instead of susceptible), which may in turn result in a different antimicrobial agent being used for treatment.

Required actions:

1. Laboratories should continue to follow their current QC procedures for Imipenem 10 µg IPM10 (ref. 065048L) for the lots listed in Table 1 in accordance with the laboratory quality assurance programme and local regulatory requirements, with a modification to increase the frequency of QC testing to every day of use.
2. Laboratories should include in the QC testing the strain *P. aeruginosa* ATCC 27853, which is defined as the stability indicator for Imipenem 10 µg IPM10 (ref. 065048L). The zone diameter for this specific strain must fall in the acceptable range to confirm the validity of the QC test and performance of the discs with the clinical isolate.
3. Only report results if all inhibition zones for QC strains are within acceptable ranges.

We apologize for any inconvenience this issue might have caused you.

We remain at your disposal for any further information.

Sincerely,
Regulatory Affairs Department

LIOFILCHEM s.r.l.
BACTERIOLOGY PRODUCTS
Via Scozia - Zona Ind.le
64026 Roseto degli Abruzzi (TE)
Cod. Fisc. e Partita IVA 00530130673

NOTE: The Competent Authorities have been informed of this communication.