Field Safety Notice

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Ponte San Nicolò, 24 Jun 2025

Field Safety Notice

Revision of the Instructions for Use of TISSUE-C (ref. MCM-246 Rev. 11)

Manufacturer: AL.CHI.MI.A. S.r.l.

Product name: TISSUE-C

Product reference: TIS 001-00

FSN Reference No: FSN000453

Dear Customer,

AL.CHI.MI.A. S.r.l. is writing to inform you of a revision to the IFU⁴ (ref. MCM-246) for the TISSUE-C medical device (product ref. TIS 001-00), a medium intended for corneal storage at 31±2°C for up to 30 days.

This revision is part of a corrective action undertaken by AL.CHI.MI.A. S.r.l. following reports received from the field in France, as previously communicated in the Field Safety Notice ref. FSN_07_2024_TIS 001-00, dated 24/07/2024.

The purpose of this notice is to inform users of the updated IFU and to emphasize the importance of following the revised instructions to ensure the safe and effective use of the TISSUE-C medical device.

The same IFU (ref. MCM-246) for TISSUE-C (product ref. TIS 001-00) is shared with CARRY-C⁵ medical device (product ref. CAR 001-00), which is intended for sequential use following TISSUE-C. For this reason, this notice also applies to the CARRY-C medical device.

Action being taken by AL.CHI.MI.A. S.r.l.

To support the correct use of the TISSUE-C medical device, **AL.CHI.MI.A. S.r.l. has implemented a revision of the product's IFU.** The revised IFU are intended to provide clear instructions on how to handle the AL.CHI.MI.A products to ensure best practices in storage of human donor corneas.

There are no changes to the design or composition of the device, which continues to meet all applicable safety and performance requirements.

A summary of the key updates included in the revised IFU is provided in **Annex 1**.

⁴ Instructions for Use

⁵ Medium for corneal deswelling and transport at room temperature for up to 5 days

Timeline for updated IFU

The revised IFU will be implemented starting with lot number P25050023 of TISSUE-C (product ref. TIS 001-00) and lot number P25030020 O.N.⁶ 250258 of CARRY-C (product ref. CAR 001-00).

Advice on action to be taken by the user

- Carefully read and understand the revised IFU for the new instructions provided.
- Ensure that all personnel involved in the handling, storage, and use of the TISSUE-C (product ref. TIS 001-00) medical device are informed of the IFU updates and trained accordingly.
- Follow the updated instructions strictly to guarantee proper use.
- Complete the attached **Customer response form (Annex 2)** and return it to AL.CHI.MI.A. S.r.l. (vigilance@alchimiasrl.com) within ten (10) days after receipt to confirm acknowledgment and understanding of this communication.

For additional information regarding this communication, please contact:
vigilance@alchimiasrl.com.

Thank you for your attention.

The National Competent Authority, ANSM, is informed of this safety action.

Yours Sincerely,

Patrick Quero

AL.CHI.MI.A. S.r.l.

Person Responsible for Regulatory Compliance Delegate

List of attachments

- **Annex 1:** Comparison between the current (ref. MCM-246 Rev. 10) and updated (ref. MCM-246 Rev. 11) IFU
- **Annex 2:** Customer response form

⁶ Order Number (of production)

Annex 1 - Comparison between the current (ref. MCM-246 Rev. 10) and updated (ref. MCM-246 Rev. 11) IFU

Changes are highlighted in red text.

IFU section	Current IFU (ref. MCM-246 Rev. 10)	Updated IFU (ref. MCM-246 Rev. 11)
INDICATIONS FOR USE	TISSUE-C, medium for corneal storage at 31°C for up to 30 days. CARRY-C, medium for corneal deswelling and transport at room temperature for up to 5 days. Medical devices are intended for use by Eye Banks and centres performing corneal excision and/or keratoplasty. Any other use must be considered improper. MEDICAL DEVICES DO NOT WARRANT THE SUITABILITY OF CORNEA FOR TRANSPLANTATION.	TISSUE-C, medium for corneal storage at 31±2°C for up to 30 days. CARRY-C, medium for corneal deswelling and transport at room temperature for up to 5 days. Medical devices are intended for use by Eye Banks and centres performing corneal excision and/or keratoplasty. Any other use must be considered improper. MEDICAL DEVICES DO NOT WARRANT THE SUITABILITY OF CORNEA FOR TRANSPLANTATION.
PACKAGING AND COMPOSITION	TISSUE-C, balanced salt solution Eagle's MEM, 2% bovine serum, sodium bicarbonate, 100 IU/ml penicillin G, 0,1 mg/ ml streptomycin, 0,25 µg/ml amphotericin B, and sodium pyruvate, pH 7,4. CARRY-C, balanced salt solution Eagle's MEM, dextran, 2% bovine serum, sodium bicarbonate, 100 IU/ml penicillin G, 0,1 mg/ml streptomycin, 0,25 µg/ml amphotericin B, and sodium pyruvate, pH 7,4. The CARRY-C package includes 12 singularly packaged sterile closures for the transport of corneas suitable for transplantation.	TISSUE-C, balanced saline solution of Eagle's MEM, 2% bovine serum, sodium bicarbonate, 100 IU/ml penicillin G, 0.1 mg/ml streptomycin, 0.25 µg/ml amphotericin B, and sodium pyruvate, pH 7,4. CARRY-C, balanced saline solution of Eagle's MEM, 6% dextran, 2% bovine serum, sodium bicarbonate, 100 IU/ml penicillin G, 0.1 mg/ml streptomycin, 0.25 µg/ml amphotericin B, and sodium pyruvate, pH 7.4. The CARRY-C package includes 12 singularly packaged sterile closures for the transport of corneas suitable for transplantation. TISSUE-C and CARRY-C devices are supplied sterile.
DIRECTIONS FOR USE	A. Corneal storage in TISSUE-C - Before use, thaw the solution in the refrigerator (2–8°C) or at room temperature. After thawing, stir the solution. - Carry out all the subsequent steps under aseptic conditions. - Make sure that the solution is at a temperature avoiding thermal shock to the cornea.	A. Corneal storage in TISSUE-C - Before use, thaw the solution in the refrigerator (2–8°C) overnight, or at room temperature (20-25°C) for up to 4 hours until complete thawing, or in a water bath at 37±2°C for up to 30 minutes. After thawing, stir the solution. An unopened bottle can be stored at 2-8°C for up to 24 hours after thawing. - Carry out all the subsequent steps under aseptic conditions. - Make sure that the solution is at a temperature avoiding thermal shock to the cornea.

	4. Open the bottle and flame the neck. Act so as to warrant the preservation of the solution sterility. 5. Rinse thoroughly the cornea with buffered saline solution and suspend it in TISSUE-C using CornealFloat. Alternatively lay it down, with the endothelium upward, on the bottom of the TISSUE-C bottle. 6. Flame the neck of the bottle and recap it. 7. At the end of the procedure, put the bottle containing the cornea in an incubator at 31°C for a maximum of 30 days.	4. Sanitize the bottle before opening it. Act so as to warrant the preservation of the solution sterility. 5. Rinse thoroughly the cornea with buffered saline solution. 6. Suspend the cornea in TISSUE-C, preferably using CornealFloat. 7. Recap the bottle. 8. At the end of the procedure, put the bottle containing the cornea in an incubator at 31±2°C for a maximum of 30 days. Note: If deswelling is necessary, place and transport the cornea in the CARRY-C medium before transplantation, following steps B and C. This will allow to restore the physiological thickness of the cornea. If cornea deswelling is not necessary, BEFORE TRANSPLANTATION, RINSE THE CORNEA THOROUGHLY WITH BUFFERED SALINE SOLUTION PSS-L⁷ FOR TWO TIMES, USING A VOLUME OF 30 ML FOR EACH RINSING, TO REMOVE ANY RESIDUES OF THE TISSUE-C SOLUTION.
	B. Corneal deswelling and transport in CARRY-C 8. Before use, thaw the solution in the refrigerator (2–8°C) or at room temperature. After thawing, stir the solution. 9. Carry out all the subsequent steps under aseptic conditions.	B. Corneal deswelling in CARRY-C 8. Before use, thaw the solution in the refrigerator (2–8°C) overnight, or at room temperature (20-25°C) for up to 4 hours until complete thawing, or in a water bath at 37±2°C for up to 30 minutes. After thawing, stir the solution. An unopened bottle can be stored at 2–8°C for up to 24 hours after thawing. 9. Carry out all the subsequent steps under aseptic conditions.

⁷ Product ref. GRS 003-00, AL.CHI.MI.A. Srl.

	10. Make sure that the solution is at a temperature avoiding thermal shock to the cornea. 11. Open the bottle and flame the neck. Act so as to warrant the preservation of the solution sterility. 12. Remove the cornea from the bottle of TISSUE-C and rinse it thoroughly with buffered saline solution. 13. Lay down the cornea gently, with the endothelium upward, on the bottom of the CARRY-C bottle. 14. Flame the neck of the bottle, recap it and store at room temperature until transport C. Transport in CARRY-C of corneas suitable for transplantation 1. Carry out all the subsequent steps under aseptic conditions. 2. Remove the cap from the bottle and discard it. 3. Cut the ring adhering to the bottle neck after cap removal and discard it. 4. Take another cap, supplied in a separate sterile pouch, and place it on the bottle by applying a pressure until you hear a «CLICK», which indicates that the cap holds the thread. 5. Screw up until reaching the endstop and the cap is no longer rotating. 6. Transport the cornea at room temperature. 7. BEFORE TRANSPLANTATION, RINSE THE CORNEA THOROUGHLY WITH BUFFERED SALINE SOLUTION FOR TWO TIMES, USING A VOLUME OF 30 ML FOR EACH RINSING, TO REMOVE ANY RESIDUES OF THE CARRY-C SOLUTION.	10. Make sure that the solution is at a temperature avoiding thermal shock to the cornea. 11. Sanitize the bottle before opening it. Act so as to warrant the preservation of the solution sterility. 12. Remove the cornea from the bottle of TISSUE-C and rinse it thoroughly with buffered saline solution. 13. Lay down the cornea gently, with the endothelium upward, on the bottom of the CARRY-C bottle. 14. Recap the bottle and store at room temperature until transport. C. Transport in CARRY-C of corneas suitable for transplantation 8. Carry out all the subsequent steps under aseptic conditions. 9. Remove the cap from the bottle and discard it. 10. Cut the ring adhering to the bottle neck after cap removal and discard it. 11. Take another cap, supplied in a separate sterile pouch, and place it on the bottle by applying a pressure until you hear a «CLICK», which indicates that the cap holds the thread. 12. Screw up until reaching the endstop and the cap is no longer rotating. 13. Transport the cornea at room temperature. 14. BEFORE TRANSPLANTATION, RINSE THE CORNEA THOROUGHLY WITH BUFFERED SALINE SOLUTION PSS-L FOR TWO TIMES, USING A VOLUME OF 30 ML FOR EACH RINSING, TO REMOVE ANY RESIDUES OF THE CARRY-C SOLUTION.
WARNINGS AND PRECAUTIONS	• Follow good medical practice. • Before use, store as indicated on labels.	• Follow good medical practice. • Before use, store as indicated on labels and instructions for use. • If a water bath is used for thawing:

	• Soon after thawing, small amounts of cryoprecipitates may appear in the solution, due to presence of serum, which will not affect the product performance – Warm and stir the solution before use to restore homogeneity. • Store the cornea at 31°C in TISSUE-C for a maximum of 30 days and at room temperature in CARRY-C for a maximum of 5 days. • Do not admix with other substances. • Single use only – Do not re-sterilize. Use of single use products that have been already used may cause cross-contamination. Re-sterilization of the product causes a change in its physicochemical and performance characteristics. • The integrity of the medical devices is guaranteed only if the package is intact – Do not use if the bottle is open or damaged. • The expiry date is printed on the product label and package – Do not use after the expiry date. • Expired solutions must be disposed of as sanitary waste according to regulations in force – Do not litter the environment. • Carry out all steps under aseptic conditions. • Do not use TISSUE-C and CARRY-C if solution is yellow, violet or cloudy.	○ do not to fully submerge or do not invert the unopened TISSUE-C and/or CARRY-C bottle, in order to minimize the risk of potential cross-contamination from the water bath; ○ do not to exceed the specified temperature (37±2°C) and time (up to 30 minutes) limits, in order to maintain the physico-chemical characteristics of the products; ○ sanitize the outside of the TISSUE-C and/or CARRY-C bottle before opening, to maintain aseptic conditions. • The specified thawing conditions apply to the intact and frozen TISSUE-C and/or CARRY-C bottle, stored at temperatures in the range of -25°C and -15°C. • Do not refreeze the TISSUE-C and CARRY-C bottle after thawing. • Soon after thawing, small amounts of cryoprecipitates may appear in the solution, due to presence of serum, which will not affect the product performance – Warm and stir the solution before use to restore homogeneity. • Store the cornea at 31±2°C in TISSUE-C for a maximum of 30 days and at room temperature in CARRY-C for a maximum of 5 days. • Do not admix with other substances. • Single use only – Do not re-sterilize. Use of single use products that have been already used may cause cross-contamination. Re-sterilization of the product causes a change in its physicochemical and performance characteristics. • The integrity of the medical devices is guaranteed only if the package is intact – Do not use if the bottle is open or damaged. • The expiry date is printed on the product label and package – Do not use after the expiry date. • Expired solutions must be disposed of as sanitary waste according to regulations in force – Do not litter the environment. • Carry out all steps under aseptic conditions. • Do not use TISSUE-C and CARRY-C if solution is yellow, violet or cloudy.
--	--	---

- During use, change of colour of the solution to violet can occur if the bottle is closed improperly, left open or open repeatedly; change of colour to yellow can occur as a result of tissue contamination. Corneal epithelial cell detachment can cause opalescence of the solution.
- The medical devices must only be used by qualified Eye Bank personnel and/or surgeons in Eye Banks and Centres performing corneal excision and/or keratoplasty. Any residues must be disposed of as sanitary waste. In case of accidental breakage or spillage before use, confine the spill, absorb it in absorbing material and dispose of in compliance with procedures in force. In case of spillage during or after use, confine the spill, immediately wash the contaminated areas with copious amounts of water and disinfectant, and sanitize the work area.

- During use, change of colour of the solution to violet can occur if the bottle is closed improperly, left open or open repeatedly; change of colour to yellow can occur as a result of tissue contamination. Corneal epithelial cell detachment can cause opalescence of the solution.
 - The medical devices must only be used by qualified Eye Bank personnel and/or surgeons in Eye Banks and Centres performing corneal excision and/or keratoplasty. **The user shall follow protocols and controls established by the Eye Banks and the national/international guidelines relating to the qualitative and microbiological controls of corneas before transplantation. For further information, please consult European Eye Banks guidelines.**
 - **Make sure that the volume of TISSUE-C during storage of the cornea is greater than or equal to 80% of the initial volume (i.e., 100 ml). If the residual volume of TISSUE-C is less than 80 ml (e.g., due to collection of samples for microbiological analysis), it is recommended to transfer the cornea into a new TISSUE-C medium.**
 - Any residues must be disposed of as sanitary waste. In case of accidental breakage or spillage before use, confine the spill, absorb it in absorbing material and dispose of in compliance with procedures in force. In case of spillage during or after use, confine the spill, immediately wash the contaminated areas with copious amounts of water and disinfectant, and sanitize the work area.
 - **For the proper storage and preservation of the cornea, the use of the CornealFloat device, included with the TISSUE-C product, is recommended.**
 - **Rinse the cornea with buffered saline solution (e.g., PSS-L). Repeat if necessary.**
- In the event of reduced visualisation of endothelial cells during microscopic analysis, apply a hypotonic solution to improve**

	• TISSUE-C and CARRY-C contain a mixture of 100 IU/ml penicillin G, 0,1 mg/ml streptomycin, and 0,25 mg/ml amphotericin B. According to the material safety data sheet, contact with the mixture can cause ocular, skin, and respiratory irritation. In addition, there is limited evidence, based upon animal studies, that the mixture may cause birth defects following ingestion. In the absence of further studies, the mixture shall be considered as a possible teratogen in humans. It is recommended to handle it with normal precautions (do not swallow, do not inhale) and to wear protecting clothing and an appropriate eye/face protection. • BEFORE POSSIBLE CLINICAL USE OF TISSUE, RINSE THE CORNEA TWICE WITH 30 ML OF BUFFERED SALINE SOLUTION FOR EACH RINSING. • TISSUE-C and CARRY-C contain 2% bovine serum. Controls on the origin of animals and tests on component warrant the absence of endotoxins, mycoplasma, viruses, and cytopathogenic agents. It is recommended to handle the solutions with normal precautions (do not swallow, do not inhale) and to wear protecting clothing and an appropriate eye/face protection. • Used medical devices are potentially bio-hazardous. Take the necessary steps for self-protection and dispose of in compliance with national regulations on workplace safety and quality assurance principles. • Dispose of empty packages in compliance with regulations in force.	visualisation of endothelial cells. In the case Trypan Blue is used for the assessment of endothelial cell density and endothelial cell mortality, it is recommended not to exceed the concentration of 0.25% and a maximum application time of 1 minute per analysis. • If deswelling is necessary to restore the physiological thickness of the cornea before transplantation (e.g., for penetrating keratoplasty), make sure to transfer it in the CARRY-C medium. • TISSUE-C and CARRY-C contain a mixture of 100 IU/ml penicillin G, 0.1 mg/ml streptomycin, and 0.25 µg/ml amphotericin B. According to the material safety data sheet, contact with the mixture can cause ocular, skin, and respiratory irritation. In addition, there is limited evidence, based on studies on humans, that streptomycin can cause birth defects if ingested. In the absence of further studies, the mixture shall be considered as a possible teratogen in humans. It is recommended to handle it with normal precautions (do not swallow, do not inhale) and to wear protecting clothing and an appropriate eye/face protection. • BEFORE POSSIBLE CLINICAL USE OF TISSUE, RINSE THE CORNEA TWICE WITH 30 ML OF BUFFERED SALINE SOLUTION PSS-L FOR EACH RINSING. • TISSUE-C and CARRY-C contain 2% bovine serum. Controls on the origin of animals and tests on component warrant the absence of endotoxins, mycoplasma, viruses, and cytopathogenic agents. It is recommended to handle the solutions with normal precautions (do not swallow, do not inhale) and to wear protecting clothing and an appropriate eye/face protection. • Used medical devices are potentially bio-hazardous. Take the necessary steps for self-protection and dispose of in compliance with national regulations on workplace safety and quality assurance principles. • Dispose of empty packages in compliance with regulations in force.
--	--	--

CONTRA-INDICATIONS	There are no known contra-indications when following occupational safety regulations and good medical practice.	There are no known contra-indications when following occupational safety regulations and good medical practice.
SIDE EFFECTS	None known.	None known.
STERILIZATION	TISSUE-C and CARRY-C are manufactured aseptically.	TISSUE-C and CARRY-C are manufactured aseptically.
STORAGE	TISSUE-C and CARRY-C , stored at temperatures in the range of -25°C and -15°C in their sealed package, are stable up to the expiry date shown on the package label. Stability studies showed that alternate cycles of freezing and thawing do not alter the physicochemical characteristics of the media.	TISSUE-C and CARRY-C , stored at temperatures in the range of -25°C and -15°C in their sealed package, are stable up to the expiry date shown on the package label.
N/A	For any report on this product, quote the lot number on package. For any comments or suggestions, go to www.alchimiasrl.com .	In case of any serious incident occurred in relation to this product, report the incident to the Manufacturer (vigilance@alchimiasrl.com) and to the competent authority of the Member State which the serious incident occurred in, quoting the lot number indicated on the label. For any further information, go to www.alchimiasrl.com .

Dear *Customer*,

please fill this report form in full block capitals according to the notice forwarded to your attention and return it to vigilance@alchimiasrl.com. Thank you.

Annex 2 - Customer response form

Subject: Revision of the Instructions for Use of TISSUE-C (ref. MCM-246 Rev. 11)

Customer details	
Organization name	
Address	
ZIP code and City	
Country	
Telephone number	
Email	

Customer action(s) undertaken	
☐	I hereby acknowledge receipt of this letter and that I have read and understood its content.
☐	The notice has been brought to the attention of all users of TISSUE-C medical device and anyone else within our organization who needs to be informed.

Date

Name (in block capitals)

Signature

Stamp