INFORMATION CLIENT

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October 16, 2025

Affected Products:

- Automated Digital Cell Morphology Analyzer DI-60

Distributor Actions:

- ☐ Inform end-users
- ☒ Install at all customer sites
- ☐ Inform Marketing and Sales Department
- ☒ Inform Service Personnel
- ☒ Inform Application Specialists
- ☐ Retrieve and distribute updated manuals
- ☐ Other:

FSCA Barcode reader issue

Background

A field Safety Notice, FSN-CV-2025-004, is submitted on the **Automated Digital Cell Morphology Analyzer DI-60**.

Affected units

Table 1. Affected units

Reference number (REF)	Serial number (SN)	Barcode Reader
CC286297 or BS737162	63220 – 70093 Or CN63220 – CN70093	All units affected
CC286297 or BS737162	60001 – 63219 Or CN60001 – CN63219	Units affected only if they have had their original barcode reader replaced with a Honeywell Vuquest.
XU-10350	K1001 – K1095	All units affected

Figure 1. Example of label with REF and SN.


Figure 2. Image of the Honeywell Vuquest barcode reader installed in system.

Description of the product problem

It was discovered that in some instances when a slide has an unreadable barcode there is a risk that the barcode reader accidentally reads the barcode of the previously processed slide, thus resulting in a misattribution of diagnostic results.

The precondition for the issue to occur is that a slide with an unreadable barcode is run after a slide with a readable barcode in a DI-60 with either a Honeywell Vuquest or Jadak FM-5 barcode reader. In this case there is a risk that the barcode reader instead reads the barcode of the slide located in the shuttle return position and thus attributes the results from the slide in fetch position to the ID of the slide in return position. If the healthcare professional performing the review is not attentive this has the potential to lead to a wrong diagnosis.

It is less likely that the problem occurs when Sysmex standard barcodes, where the barcode is smaller and located to the left, are used.

Action to mitigate the risk by the user

1. If your unit is affected, please implement the following interim measure to mitigate the risk of sample misattribution:
 - a. Before signing each order, confirm that the number of slides in the order matches the number of slides ordered. No order shall contain more slides than prescribed.
 - b. Await contact from Sysmex service personnel, who will arrange installation of the updated CCU firmware.
2. You may continue to use the system in accordance with its intended use by following the above recommendation.

Action to mitigate the risk by service personnel

A new Firmware CCU4 version 4.2.7 has been developed to solve the problem associated with the barcode reader. This new version shall be installed in affected units through a Field Safety Corrective Action. The installation will be performed by Sysmex service personnel, who will contact affected users to schedule a time for firmware update.

Corrective actions

A new firmware for the CCU4 have been released.

1. Download the CCU Firmware for CCU4 version 4.2.7 from the CellaVision home page
2. Unzip downloaded file
3. Run the **CcuUpgradeApp.exe** to upgrade the firmware
4. Restart the system after upgrade
5. Perform **Barcode reader positioning** in the Service Wizard

Affected systems see Table 1.

Systems that have changed to a spare part barcode reader, Honeywell Vuquest, Figure 2 are affected. Honeywell Vuquest have been delivered as spare part from May 2022.

Systems with SN **60001 – 60128**. These systems were delivered with a CCU3 and if they have had their original barcode reader replaced with a Honeywell Vuquest, the CCU3 need to be replaced with CCU4 free of charge. Please contact service@cellavision.com for further assistance.

Affected DI-60s and barcode readers in stock**DI-60s**

For affected DI-60s in stock. These may be installed, but service personnel performing the installation must also upgrade to firmware version 4.2.7.

Spare part Honeywell Vuquest barcode reader

For spare part Honeywell Vuquest barcode reader, 100 103 011, in stock. These may be installed on DI-60, but Service personnel performing the installation must also upgrade to firmware version 4.2.7.

The spare part 100 103 011 Barcode reader from batch 2509291 will have an updated instruction which requires CCU4 firmware version 4.2.7.

Spare part CCU4 in stock

For spare part CCU4, 100 106 022, in stock. These may be installed, but service personnel performing the installation must also upgrade to firmware version 4.2.7. From CCU4 serial number H06320, the spare part CCU4 comes with firmware version 4.2.7 pre-installed.

For more information, contact service@cellavision.com