Abbott Diagnostics Korea, Inc.
65, Borahagal-ro, Giheung-gu,
Yongin-si, Gyeonggi-do,
Republic of Korea (16878)
T: +82.31.899.2800 / 899.2805
fax : +82.31.899.2842

Doc. No. ADK-Q2025-286

November 6, 2025

Subject: User Guidance on False Negative Results for Bioline™ TB Ag MPT64

Manufacturer: Abbott Diagnostics Korea Inc.

Country of Occurrence: France

To Whom It May Concern,

We, Abbott Diagnostics Korea Inc. (hereafter referred to as “ADK”), would like to inform you of important considerations regarding the interpretation of results obtained using the Bioline™ TB Ag MPT64. This communication is issued at the request of the ANSM following an increase in the number of false negative results related to this device reported to the ANSM. No increase in false negative reports has been observed outside of France.

Given the importance of accurate tuberculosis diagnosis, especially in high-burden or surveillance settings, it is essential to understand the potential causes of false negative results and to strictly follow the instructions for use (IFU) provided with the product.

Potential Causes of False Negative Results

1. Low MPT64 Antigen Concentration

If the MPT64 antigen concentration is below the detection limit of the product Bioline™ TB Ag MPT64, false negative results may occur.

2. Operational or Procedural Errors

Deviations from the recommended protocol may result in false negative results. Examples include:

- Delayed testing after opening the foil pouch
- Insufficient colony suspension (< 3–4 colonies in 200 µL extraction buffer)
- Incorrect sample volume (≠ 100 µL)

- Reading results outside the recommended 15 minutes window
- Misinterpreting faint test lines as negative
- Exposure of the kit to temperatures outside the recommended range (1 – 30 °C)

Recommendations to be followed by Users

Users are recommended to remain particularly vigilant regarding possible false negative results, particularly when a negative result occurs in the presence of a strong suspicion of tuberculosis. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

In this case, it is recommended to implement the following measures:

- Consider other clinically available tests if questionable results are obtained.
- Strictly follow the instructions for use provided in the Bioline™ TB Ag MPT64.
- Interpret faint test lines as positive to avoid misdiagnosis.
- Retain this letter for your laboratory records and quality assurance documentation.

ADK remains committed to providing high-quality diagnostic products and support services to meet the needs of your laboratory, healthcare providers, and patients. Should you or any of your colleagues have questions regarding this information, please contact your local Customer Service representative. We welcome further dialogue and collaboration.

Sincerely,

Signed by Taekjin Lee
 Taekjin Lee | I approve this document
11/6/2025 | 17:38 KST
ED7F79E83AE045C1891C67C971848BD9
Taekjin, Lee

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