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İMHİ

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Exmoor Plastics Ltd, 1 Western Avenue, Matrix Park,
Chorley, PR7 7NB, United Kingdom

T +44 (0) 1204 529494
E ASKVCi@Vernacare.com
www.exmoorltd.com

Advena Ltd. Tower Business Centre, 2nd Flr.,
Tower Street, Swata, BKR 4013 Malta




AURAL VENTILATION TUBES USER INSTRUCTIONS

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USER INSTRUCTIONS AURAL VENTILATION TUBES

This product may only be used, or sold, by, or on the order of, a physician

This 'grommet', otherwise known as a 'ventilation tube' or 'PE tube', is a single use device and has been sterilised by steam and is ready for use.

Do not reuse in whole or in part. Any reuse will be deemed misuse and may result in infection of a patient, user or third party. Do not use if "use by" date is expired. Refer to Unit/Lot Number label. Do not use if sterile pack is breached, damaged or contaminated.

INDICATIONS FOR USE

In the treatment of:


1. otitis media with effusion (glue ear)
2. eustachian tube dysfunction
3. other conditions requiring ventilation of the middle ear cleft, e.g. retraction pockets in the tympanic membrane, adhesive otitis

TECHNIQUE

An incision is made in the tympanum. The grommet is grasped in forceps or other insertion instrument, and the inner flange is passed through the incision into the middle ear cleft. The grommet is then manipulated so that the narrow part of the tube rests within the incision.

Identification of the inner flange

In the following illustrations, the inner flange is indicated. In those cases where the flanges are identical, either may be perceived as the inner flange, of course. If there is only one flange it is always the inner flange.

INTENDED

PERFORMANCE

This device is intended to ventilate the middle ear cleft.

RISKS/SIDE EFFECTS

Clinical

There is a small risk of infection and/or tympanosclerosis. There exists also a very slight risk that the device will be rejected or that there will be a residual perforation. There is a risk that certain drugs may increase the risk of perforation. The device may require surgical removal. There is a risk that the device may migrate into the middle ear. There is a risk that more than one Aural Ventilation Tube procedure may be required.

USER INSTRUCTIONS AURAL VENTILATION TUBES

RISKS/SIDE EFFECTS

Clinical (continued)

There is a risk of patient developing a discharge after the procedure. There is a higher risk of complications for patients with, for instance, a cleft palate. There is a risk of device being reused and / or resterilised for reuse.

General

None

HANDLING

Pack contents sterile if packaging unbreached. Do not use if packaging breached. Do not resterilise. Open the packaging with great care, in order to avoid spilling the device inadvertently, and remove the grommet aseptically.

STORAGE

Store in a dry environment. Avoid extremes of temperature and humidity. Optimal storage conditions are between 10 - 35°C and 20 - 80% relative humidity. Keep out of direct sunlight.

PATIENT BRIEFING

Ears should be kept dry for the first four weeks. The patient should be advised to submit to regular supervision by his/her physician until the tube has extruded, or been removed, from the ear.

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