Date: 2025-12-03

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For Attention of*: Vigilance/Quality manager of the facility products.

Dear customer,

We would like to inform you about a precautionary measure concerning the product Soluscope™ Serie 5 washer-disinfector with the following serial numbers: 30414, 30415, 30233, 30234, 30245, 29243, 29244, 29245, 30248, 30417 and 29239.

As part of ongoing quality checks, we have identified a slight gap between the connector block and the bowl wall in some Soluscope™ Serie 5 units. This gap could potentially allow contaminants to accumulate and later be released into endoscopes during reprocessing. To date no incident has been reported. As a precautionary measure, we are issuing the attached Field Safety Notice to inform customers and advise them that a technician will visit shortly to verify the situation and initiate the necessary corrective actions.

We kindly ask you to take the following actions:

- Please do not use or distribute the identified products. If you have distributed these products to your customers, it is essential that you inform them about the issue and advise them to stop use or distribution.
- Complete the attached reply form: review the information in this document, follow the appropriate actions outlined in section 3, and confirm your understanding and compliance with the outlined measures. Please return the Customer Reply Form to us within four weeks after receiving the FSN. In case you have distributed the products to your customers, please collect their responses and complete a single reply form with aggregated data.

We sincerely apologize for any inconvenience this may cause and appreciate your understanding and cooperation in this matter.

Thank you for your cooperation and understanding.

Best regards,
ECOLAB VIGILANCE
On behalf of Soluscope SAS

Field Safety Notice (FSN)

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Non-compliant bowl

1. Information on Affected Devices*	
1.	1. Device Type(s)*
	Class IIb medical device under MDR
1.	2. Commercial name(s)*
	Soluscope™ Série 5
1.	3. Unique Device Identifier(s) (UDI-DI)
	359761SS2B54S5EQPT
1.	4. Primary clinical purpose of device(s)*
	The Soluscope™ Série 5 automated washer-disinfector (WD) is intended to clean and disinfect semi-critical, heat sensitive flexible devices with or without channels. It is intended for use exclusively with its dedicated, single-use cleaner Soluscope ECL-Enzymatic Cleaner and disinfectant Soluscope PAA. It is destined to be used by trained personnel familiar with device reprocessing, handling cleaner(s) and disinfectant(s), in a hospital or medical setting, in endoscopy departments, operating theaters or medical offices.
1.	5. Device Model/Catalogue/part number(s)*
	Soluscope™ Serie 5 - SL S5
1.	6. Software version
	N/A
1.	7. Affected serial or lot number range
	30414, 30415, 30233, 30234, 30245, 29243, 29244, 29245, 30248, 30417 and 29239
1.	8. Associated devices
	N/A

2. Reason for Field Safety Corrective Action (FSCA)*	
2.	1. Description of the product problem*
	As part of ongoing quality checks, we have identified a slight gap between the connector block and the bowl wall in some Soluscope™ Serie 5 units.
2.	2. Hazard giving rise to the FSCA*
	To date, no incident has been reported. However, this gap could potentially allow contaminants to accumulate and later be released into endoscopes during reprocessing.

3. Type of Action to mitigate the risk*		
3.	1. Action To Be Taken by the USER* <div style="display: flex; justify-content: space-between;"> <div> <input checked="" type="checkbox"/> Identify Device <input type="checkbox"/> Return Device <input checked="" type="checkbox"/> On-site device modification / inspection <input checked="" type="checkbox"/> Inform impacted users within your facility <input type="checkbox"/> Other </div> <div> <input type="checkbox"/> Quarantine Device <input type="checkbox"/> Destroy Device <input type="checkbox"/> None </div> </div> <p>Provide further details of the action(s) identified.</p>	
	2. Action To Be Taken by the DISTRIBUTOR* <div style="display: flex; justify-content: space-between;"> <div> <input checked="" type="checkbox"/> Identify Device <input checked="" type="checkbox"/> Return Device OR <input checked="" type="checkbox"/> On-site device modification / inspection <input checked="" type="checkbox"/> Inform all your impacted customers <input type="checkbox"/> Other </div> <div> <input type="checkbox"/> Quarantine Device <input type="checkbox"/> Destroy Device <input type="checkbox"/> None </div> </div> <p>Provide further details of the action(s) identified.</p>	
3.	3. By when should the action be completed?	Immediately
3.	4. Is customer Reply Required? * (If yes, form attached specifying deadline for return)	Yes Within 4 weeks after receipt of this FSN
3.	5. Action Being Taken by the Manufacturer* <div style="display: flex; justify-content: space-between;"> <div> <input checked="" type="checkbox"/> Product Removal <input type="checkbox"/> Software upgrade <input type="checkbox"/> Other </div> <div> OR </div> <div> <input checked="" type="checkbox"/> On-site device modification/inspection <input type="checkbox"/> IFU or labelling change <input type="checkbox"/> None </div> </div>	

4. General Information*		
4.	1. FSN Type*	New
4.	2. Manufacturer information (For contact details of local representative refer to page 1 of this FSN)	
	a. Company Name	Soluscope SAS
	b. Address	100, rue Fauge – Z.I. les Paluds, 13400 Aubagne – France
	c. Website address	www.soluscope.com
4.	3. The Competent (Regulatory) Authority of your country has been informed about this communication to customers. *	
4.	4. List of attachments/appendices:	Annex A: Customer reply form

Transmission of this Field Safety Notice	
	<p>This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)</p> <p>Please transfer this notice to other organisations on which this action has an impact. (As appropriate)</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback.</p>

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