Boston Scientific

Boston Scientific Corporation
4100 Hamline Ave North
St. Paul, MN 55112-5798 USA
+1.651.582.4000
SRN US-MF-000010993
www.bostonscientific.com

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URGENT Product Advisory - Field Safety Notice (FSN) Reference: 97125289H-4-FA

Information on Affected Devices

This information applies to all listed models. No serial number identification required.

Models: L121 / L131 (ESSENTIO™), L221 / L231 (PROONENT™), L321 / L331 (ACCOLADE™), S722 (ALTRUA™), U125 / U128 (VALITUDE™), U225 / U226 / U228 (VISIONIST™)

Reason for Field Safety Corrective Action

Any device that experiences wandless ZIP™ telemetry disablement by the battery impedance test due to detection of high battery impedance requires replacement before normal battery replacement is indicated.

There is a 7.6% likelihood that an individual CRT-P or DR-EL will need to be replaced early due to high battery impedance induced ZIP disablement. For those devices, the projected reduction in longevity is $10.9\% \pm 9.6\%$ (mean \pm standard deviation). 92.4% of CRT-P and DR-EL devices are expected to achieve anticipated longevity; thus, the overall longevity impact is 1%.¹

Recommendations**Consider these updated longevity expectations when making implant decisions.**

Please forward this notice to all relevant personnel within your organization.

The Regulatory Authority of your country has been informed of this Field Safety Notice. Adverse events should be reported to Boston Scientific and to your national Regulatory Authority, as applicable.

¹ Longevity reduction estimates are based on data included in the March 2026 Brady SMR6 letter. Any updates to these data will be published in the Product Performance report at www.bostonscientific.com/ppr.

Additional Information

Up-to-date product performance information, including this topic, and a device lookup tool is available within our Product Performance Resource Center at www.bostonscientific.com/ppr.

If you have additional questions regarding this information, please contact your Boston Scientific representative or Technical Services.

Alexandra Naughton
Vice President, Quality Assurance

United States Technical Services

1.800.CARDIAC (227.3422)
tech.services@bsci.com

International Technical Services

+32 2 416 7222
intltechservice@bsci.com

Asia Pacific Technical Services

+61 2 8063 8299
aptechservice@bsci.com

The following are trademarks of Boston Scientific or its affiliates: ACCOLADE, ALTRUA, ESSENTIO, PROPONENT, VALITUDE, VISIONIST, ZIP. All other trademarks are the property of their respective owners.

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All Rights Reserved.
© 2026 Boston Scientific Corporation or its affiliates.