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To the attention of laboratory managers, local reactivigilance correspondents and directors of the health establishment

SUBJECT: SAFETY INFORMATION

ANSM file reference: R2429421

Device Type: MTS™ Temocillin TMO 0.064-1024

Ref.: 920291, 92029, 920290

Lot: 032823060, 021924042

Roseto degli Abruzzi, 10 December 2024

Dear Customer,

This letter is intended for users of the following MTS product:

Temocillin TMO 0.064-1024 (ref. 920291, 92029 and 900290).

Description of the issue

Several customers have reported out-of-range high MIC results for *Escherichia coli* ATCC 25922 when performing QC testing with Temocillin TMO 0.064-1024 (lot 032823060 and 021924042).

Liofilchem internal investigation identified that:

- QC tests on the retained samples of the impacted lots of Temocillin TMO 0.064-1024 complied with specifications.
- QC failures were observed depending on the Mueller Hinton medium supplier used.
- There is a potential performance issue for QC strain (*E. coli* ATCC 25922) and clinical isolates on Temocillin TMO 0.064-1024 (whatever the lot) used in association with cation-adjusted Mueller Hinton Agar (MHA) other than that from Liofilchem and BBL/BD that could lead to False Resistant results. A note regarding the agar medium is also included in the actual batch related certificate of analysis (available from liofilchem.com/coa).

Clinical impact

No clinical impact for the patient. As per our IFU, quality control test must be performed to check the performance of MTS and patient isolate results should not be reported if the QC results are outside of specifications. In case performance characteristics are not checked, the only risk is that there is a potential to discard the molecule due to erroneous categorisation (resistant instead of susceptible) as other therapeutic alternatives are possible.

Required actions:

1. Perform systematically the QC with *E. coli* ATCC 25922 to detect performance issues.
2. If QC fails (MIC out of range high), repeat the test using a different brand/lot of MHA plates than currently used.
3. When results with another MHA brand are the same, contact immediately Liofilchem or the local representative for assistance.

We apologise for any inconvenience this issue might have caused you.

We remain at your disposal for any further information.

Sincerely,
Regulatory Affairs Department

LIOFILCHEM s.r.l.
BACTERIOLOGY PRODUCTS
Via Scozia - Zona Ind.le
64026 Roseto degli Abruzzi (TE)
Cod. Fisc. e Partita IVA 00530130673

NOTE: The ANSM has been informed of this communication.