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7 septembre 2026

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
Date/Signature : _____

Tracoe Vario

**Specific Attention to Details
When Using Tracoe Vario Tracheostomy Tubes
with Minimally Traumatic Insertion System**

REF 450-P, 451-P, 460-P, 461-P, 470-P, 471-P

Contents	
EN - ENGLISH.....	4

Image A



Image B

EN - ENGLISH

Specific Attention to Details When Using Tracoe Vario Tracheostomy Tubes with Minimally Traumatic Insertion System

REF 450-P, 451-P, 460-P, 461-P, 470-P, 471-P

Important Note: This supplement does not replace the original Instructions for Use (IFU). Read the applicable IFU carefully before use. The IFU forms an integral part of the product documentation and must be available at all times. For the safety of both patients and users, pay particular attention to the following information.

Specific Attention When Preparing the Minimally Traumatic Insertion System

CAUTION

- Appropriate measures must be available to ensure immediate ventilation through an alternative airway in the event of complications during the procedures described (e.g. translaryngeal intubation or laryngeal mask airway).
- **Do not disassemble the pre-assembled tracheostomy tube and insertion system** consisting of the tracheostomy tube, inserter and guiding catheter.
- **Do not pull on the tip of the guiding catheter or silicone sleeve solely.**

1. Visually inspect the device and verify that the insertion system is correctly assembled as described in the applicable IFU.
2. Lubricate the silicone sleeve and the distal end of the tracheostomy tube using a pea-sized amount of the supplied lubricating jelly.

Specific Attention When Using the Tube with the Minimally Traumatic Insertion System

CAUTION

- **Do not pull on the guiding catheter separately while the inserter is still inside the tube.**
- **Do not bend the tracheostomy tube during the removal of inserter, catheter and Seldinger wire.**

After tube insertion, withdraw the orange inserter first, followed by the guiding catheter and Seldinger wire.

1. Stabilize the tracheostomy tube with one hand. Withdraw the inserter only, leaving the guiding catheter and Seldinger wire in place (**see page 4, Image A**).
2. After removing the inserter, withdraw the guiding catheter together with the Seldinger wire while maintaining the tracheostomy tube in its position within the trachea (**see page 4, Image B**).

- After removing the guiding catheter, ensure that the silicone sleeve is still located on the guiding catheter. If the silicone sleeve detaches from the guide catheter, proceed as follows:

- Immediately inspect the lumen of the tracheostomy tube. If the sleeve remains in the tube, it must be replaced with the spare device immediately.
- If the silicone sleeve cannot be visually identified within the lumen of the tube, obtain an immediate chest X-ray to determine its location. As the silicone sleeve is radiopaque, radiographic imaging facilitates subsequent bronchoscopic retrieval.

2026-09
801590-03

Tracoe Medical GmbH
Reichelsheimer Straße 1 / 3
55268 Nieder-Olm, Germany
info-tracoe@atosmedical.com
www.atosmedical.com

Atos and the Atos Medical logo
are trademarks of Coloplast A/S.
© 2026 Coloplast A/S.
All rights reserved.